Leading in Integrated Beverage Solutions

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*Top quality drinks require consistently the highest product safety standards and this necessitates a constant review and monitoring regime. Döhler **Micro-safety Design**[®] (DMD) guarantees this standard and provides a platform for excellent food safety practice along the entire added-value chain. Döhler, in collaboration with Professor Werner Back of Munich/Weihenstephan Technical University, has been leading the nutrient-media field of microbiological monitoring for the soft drinks and brewing industries for many years.*

*Using the Döhler detection media for example NBB[®], OFS-Agar, TransFast, LMC etc. and modern molecular biological methods (*Vermicon Identification Technology*, VIT) for simple, quick and efficient detection of microbiological contaminations for all beverage types,*

microorganisms can be identified long before sub-standard finished goods leave the factory premises. Virtually all Döhler detection media are provided ready-to-use.

Meanwhile, the Bypass-Membrane-System (BM-System) is a special sampling technique designed for production plant validation. For special tailor-made customer analysis Döhler is also able to offer a bespoke service based in their own laboratories.

Our long-time experience together with the special knowledge of our partners in the Döhler “Innovation Network” provides pragmatic solutions for microbiological questions – flexible, rapid and reliable.

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Detection of Microorganisms: Reliable and Fast

The Döhler nutrient media for the soft drinks industry give the user the possibility to detect harmful microorganisms reliable and fast. The culture media are very selective: only spoiling microorganisms will have optimal conditions for growth. All culture media are ready-to-use, additional mixing and sterilization are not necessary.

The inspector's stamp of the Technical University of Munich/Weihenstephan guarantees controlled quality.

*With OFS-Agar, SSL-Broth and the TransFast-System (Broth and Gel) harmful moulds, yeast, lactic and acetic acid bacteria could be verified. TransFast is characterized by a special rapidness. For pre-enrichment of harmful microorganisms TransFast-Broth will be the best choice. Specially for analyzing water LMC-Concentrate and DEV-Nähragar are available for detection of *E. coli* and coliform bacteria respectively Total Viable Count (TVC).*

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Nutrient Media for the Isolation of Spoiling Microorganisms in Beverages

(Yeasts, moulds, acetic and lactic acid bacteria)

Application:

- Soft drinks
- Dairy
- Winery

- a) *Pour plate method, Ready-to-Drink (RTD) beverages, rinse water samples, cultivation of microorganisms*
- b) *Nutrient base for membrane filters*

Preparation:

(pour plate method)
Place approx. 3 ml or 10–15 ml of the test sample into the Petri dishes (9 cm diameter and 14 cm diameter respectively) and fill with warm agar to a height of

3–4 mm. Swirl the mixture in the Petri dishes. Place test sample into an incubator for three days at a temperature of 28 °C.

Analysis:

A good evaluation of the results can be achieved by the macroscopic examination of samples due to the individual growth characteristics and odorous emissions of each species.

Characteristics:

Yeasts

Mat, yellow-whitish or red colonies. Partially star-like growth or formation of a white film (non-fermenting yeast).

Moulds

Form large, circular mycel with white edge (diameter usually above 0.5 cm) with green, brown or black centre.

Lactic acid bacteria

Small, shiny, elliptical colonies on the surface of the nutrient media. Odour: slightly acidic, resemble rancid butter. Catalase-negative.

Acetic acid bacteria

Growth limited to surface. Shiny; extremely pungent acetic acid odour. Catalase-positive.

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Isolation of Spoiling Microorganisms in Beverages

(Yeasts, moulds, acetic and lactic acid bacteria)

Application:

- Soft drinks
- Dairy
- Winery

Tracer studies in filterable and non-filterable beverages and base ingredients.

Test medium of swabs, cultivation of spoiling microorganisms of RTD beverages.


Preparation:

Mix 200 ml of RTD beverages resp. 10–20 ml concentrate, puree, compound or sugar syrup with 50 ml broth and keep in an incubator for 2 days at 28 °C.

Prepare pour plate with 3 ml test medium (culture medium OFS-agar). For this it is sufficient to incubate the sample for 2 days at 28 °C as the microorganisms are at their optimum growth phase due to enrichment.

Analysis:

The interpretation of the results is primarily a qualitative one. With some experience it is, however, possible to estimate the initial number of organisms that resulted in the final number of microorganisms (see diagram).

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Nutrient Agar as per DEV (pH 7.3 ± 0.2)

Application:

Evaluation of the Total Viable Count (TVC) in water as per German standard method (DEV) and drinking water regulations.

Evaluation of Total Viable Count (TVC) in various food products and rinse water samples. Pour plate method, nutrient base for membrane filters.

Preparation:

For water analysis, either prepare 1 ml of the test sample according to the pour plate technique or filter the appropriate quantity (using a membrane filter).

Water analysis as per DEV: incubate 1 water sample for 24 hours at 37 °C, incubate 1 water sample for 48 hours at 20 °C (maximum permitted number of microorganisms: 100 cfu per ml at 20 °C; 20 cfu per ml at 37 °C).

For evaluation of total number of microorganisms from different test materials the appropriate actual regulations apply.

Analysis:

After an appropriate period of incubation (usually 2 days) the colonies are counted. If the number of microorganisms can not be established testing must be repeated using a serial dilution of the sample.

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DEV-Lactose- Broth (pH 7.2 ± 0.1)

Application:


Qualitative presumptive test for evaluation of *E. coli* and coliform bacteria in water.

Preparation:

Add test sample (250 ml) to concentrate in test bottle and incubate for a maximum of 48 hours at 37 °C (see graphic). Using the method for carbonated mineral water, the water must be decarbonated before testing (e. g. decarbonation by shaking or neutralisation with sodium hydroxide).

Analysis:

If there is a definite change in colour to (opaque) yellow, cloudiness or gas is released by the sample during incubation, the result is positive. A more definite identification is possible by addition to Endo-agar or similar nutrient media and subsequent colour sequence test.

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Handling of the Nutrient Media OFS-Agar and DEV Nutrient Agar

The unopened bottle containing solid nutrient media will be heated at max. 95 °C until the media is liquefied completely. Let the media cool down slowly to a temperature of 45 to 48 °C. Differences of temperature more than 20 °C must be avoided.

Swab Samples:

As test medium for sterile swab applications in test tubes SSL-broth or physiological saline solution can be used. Swabs taken from the test medium are used for a pour plate with appropriate agar.

Catalase-Test:

This is a test to differentiate between acetic acid and lactic acid bacteria. Apply one drop of a 5 % hydrogen peroxide solution (H₂O₂) to a single colony. The Catalase-positive acetic acid bacteria produce a strong effervescent reaction, whereas the lactic acid bacteria show no reaction at all.

GRAM-Screening (KOH-Test):

Instead of the complex GRAM-staining-test it is possible to evaluate bacteria with the so-called "KOH-test". This screening is based on different cell wall structures of GRAM-positive and GRAM-negative microorganisms.

A suspension of microorganisms (pure culture) is treated with caustic potash solution (5 %). Preferably take one big colony or collect a number of small, similar colonies from Petri dishes and

mix these material on a microscope slide by an inoculation loop together with a drop of caustic potash solution. After a short time GRAM-negative bacteria will show slimy threads if dabbed with an inoculation loop. GRAM-positive bacteria will show no reaction.

Oxidase Detection:

The detection of the enzyme Oxidase is one of the key tests for differentiation of Enterobacteria (Oxidase-negative) and other non-fermenting GRAM-positive bacteria, which are normally Oxidase-positive.

For evaluation of the colonies in a fresh culture (Petri dish) a few drops of 1 % aqueous solution of N,N-Dimethyl-1,4-phenylenediamid-dihydrochlorine are added. Preferably, the nutrient media should be free of glucose and nitrate. Colonies of Oxidase-positive microorganisms turn pink to crimson, later dark brown to black. Ready-to-use test kits are also available. In this case a suspension of microorganisms will develop colours on a test stripe in the same way as described above.

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TransFast – Rapid and Selective Detection of Microorganisms Harmful to Soft Drinks

This test is based on a broth which concentrates non-filterable ready-to-drink beverages prior to examination for



microbe content and enables a larger volume (approx. 200 ml) to be analyzed. The second culture medium is a gel which enables membrane filters, rinse water samples, ready-to-drink beverages, raw materials and the concentrated beverages mentioned above to be examined directly.

The test detects microorganisms which can grow in a pH-environment of less than 4.5 (beverages), and spoil the product (mould, yeast, lactic and acetic acid bacteria). Microbiological results can usually be obtained after 24 hours or at most two days. An additional pre-enrichment procedure requires again 1–2 days.

Preparation:

Place the material to be examined in a sterile 80 ml test tube and fill with appropriate amount of gel under sterile conditions.

Place the samples in a glass-fronted viewing box which can be kept at an incubating temperature of 28 °C.

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
Finally turn the bottle upside down to ensure that the headspace is sterilized. The prepared culture medium is pale opal in colour and may be streaky. However, this will not effect any tests.

Shelf life of the prepared culture medium

After cooling, the culture medium can be kept for up to 4 weeks at 5 °C. After 4 weeks it should be used up or discarded.


TransFast – Making up the Culture Media

pH-adjustment

When the solution reaches 40–50 °C adjust the pH-level to 4.3–4.4 using component 2. Use a pH-meter which compensates for temperature to carry out the adjustment.

Pasteurizing the culture medium

When the pH-level has been adjusted, heat the solution to 80 °C, and maintain it at this temperature for approximately 60 seconds. Pour it into prepared, clean bottles and seal. The bottles and screw tops must be able to withstand a temperature of 80 °C.

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TransFast – Analysis

After one or two days, any cultures in the gel can be detected with the naked eye. Even at this stage, the cultures which have grown can be distinguished macroscopically by an experienced technician.

However, as in the Petri dish method, the colonies must be examined under a microscope to identify them more accurately.


The relatively low pH-value in the culture medium inhibits growth and only allows a selective number of cultures, which cause spoilage, to grow.

It is therefore quite possible that certain cultures will not grow in the gel, particularly in rinse water samples, which can otherwise be detected in the classic Petri dish method.

However, they have absolutely no spoilage effect on the ready-to-drink beverage. They can be ignored, and have no bearing on the analysis of microbes found in the gel. Equipment is generally washed with normal water which has not been pasteurized, and consequently the microbes found in water are present. These cannot be detected and do not have any effect.

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Test sample	Location of test points	Frequency of test samples	Microorganisms to be examined	Nutrient media			
				DEV lact. broth	DEV stand- ard	OFS TFG	SSL TFB
1 Water	After well pump	2 x per month	Total viable count, coliform bacteria	•	•		
2 Water	After mains water flow meter	2 x per month	Total viable count, coliform bacteria	•	•		
3 Water	After pressure tank or reservoir tank	2 x per month	Total viable count, coliform bacteria	•	•		
4 Water	After water treatment and water flow meter	2 x per month	Total viable count, coliform bacteria and beverage spoiling organisms	•	•	•	
5 Water	In syrup chamber before mixing vessel	2 x per month	Beverage spoiling organisms			•	
6 Water	Mixer inlet	1 x per week	Beverage spoiling organisms			•	•
		1 x per month	Total viable count, coliform bacteria	•	•		
7 Water	Watermixer outlet (for mineral water)	1 x per week	Total viable count, coliform bacteria	•	•		
8 Water	Filling inlet	1 x per week	Total viable count, coliform bacteria	•	•		
9 Water	From filled bottles	2 x per bottling operation	Total viable count, coliform bacteria	•	•		
10 Sugar syrup	From tanker	each delivery	Beverage spoiling organisms			•	
11 Rinse water	From storage tank	before each bottling operation	Beverage spoiling organisms			•	
12 Base material	From packaging	each delivery	Beverage spoiling organisms			•	
13 Base material	From inlet hose to mixing vessel	2 x per week	Beverage spoiling organisms			•	
14 Final syrup starting solution	From mixing vessel	2 x per week	Beverage spoiling organisms			•	•
15 Final syrup	Before mixer	2 x per week	Beverage spoiling organisms			•	•
16 Final beverage	Mixer outlet	2 x per week	Beverage spoiling organisms			•	•
17 Final beverage	From filled bottles	2 x per filling operation	Beverage spoiling organisms			•	•
18 Rinse water	Mixing vessel	after cleaning	Beverage spoiling organisms			•	
19 Rinse water	Inlet to mixer	after cleaning	Beverage spoiling organisms			•	
20 Rinse water	Mixer outlet	after cleaning	Beverage spoiling organisms			•	
21 Rinse water	From filler	after cleaning	Beverage spoiling organisms			•	
22 Cleaned bottle	After rinser	1 x per week 5 bottles	Total viable count, coliform bacteria and beverage spoiling organisms	•	•	•	
23 Post-spray cold water	Rinser	1 x per week	Total viable count, coliform bacteria and beverage spoiling organisms	•	•	•	
24 Post-spray hot water	Rinser	1 x per week	Total viable count, coliform bacteria and beverage spoiling organisms	•	•	•	
25 Swab samples	Bottle outlet at rinser	1 x per week	Total viable count, coliform bacteria and beverage spoiling organisms	•	•	•	
26 Swab samples	Bottle inspection	1 x per week	Total viable count, coliform bacteria and beverage spoiling organisms	•	•	•	
27 Swab samples	Filler valves	1 x per week	Total viable count, coliform bacteria and beverage spoiling organisms	•	•	•	
28 Swab samples	Filler: broken pieces spray	1 x per week	Total viable count, coliform bacteria and beverage spoiling organisms	•	•	•	
29 Swab samples	Filler casing	1 x per week	Total viable count, coliform bacteria and beverage spoiling organisms	•	•	•	
30 Swab samples	Capper: point of capping	1 x per week	Total viable count, coliform bacteria and beverage spoiling organisms	•	•	•	

Process Control One-Way Bottle

Test sample	Location of test points	Frequency of test samples	Microorganisms to be examined	Nutrient media			
				DEV lact. broth	DEV stand- ard	OFS TFG	SSL TFB
22 Rinsed bottles	After rinser	1 x per week 5 bottles	Total viable count, coliform bacteria and beverage spoiling organisms	•	•	•	
23 Rinse water	Before rinser	1 x per week	Total viable count, coliform bacteria and beverage spoiling organisms	•	•	•	
24 Swab samples	Water discharge rinser	1 x per week	Total viable count, coliform bacteria and beverage spoiling organisms	•	•	•	
25 Swab samples	Gripper	1 x per week	Total viable count, coliform bacteria and beverage spoiling organisms	•	•	•	

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Microbiological Culture Media for Practical Operational Monitoring

Product	Article-No.	Package
OFS-Agar Application: <ul style="list-style-type: none"> • Soft drinks, dairy, winery • Ready-to-drink (RTD) beverages, rinse water samples, cultivation of microorganism • Pour plate method, nutrient base for membrane filters 	204707782	9 x 250 ml
SSL-Broth Application: <ul style="list-style-type: none"> • Soft drinks, dairy, winery • Tracer studies in filterable and non-filterable beverages and base ingredients • Test medium of swabs, cultivation of spoiling microorganism of RTD beverages 	204712782	9 x 250 ml
LMC-Concentrate Application: <ul style="list-style-type: none"> • Qualitative presumptive test for evaluation of E.coli and coliform bacteria in water 	204713700	9 x 50 ml
DEV-Agar Application: <ul style="list-style-type: none"> • Evaluation of the Total Viable Count (TVC) in water as per DEV and drinking water regulations • Evaluation of the Total Viable Count (TVC) in various food products and rinse water samples • Pour plate method, nutrient base for membrane filters 	204726782	9 x 250 ml
TransFast-Broth (TF-B) Application: <ul style="list-style-type: none"> • Broth which concentrates non-filterable ready-to-drink beverages prior to examination for microbe content and enables a larger volume (approx. 200 ml) to be analyzed 	204727782	9 x 250 ml
TransFast-Gel 1 (TF-G 1 Powder) Preparation: TF-G 1 + TF-G 2 Application: <ul style="list-style-type: none"> • Gel which enables membrane filter, rinse water samples, ready-to-drink beverages, bases, raw materials and the concentrated beverages mentioned above to be examined directly • Microorganism which can grow in a pH-environment of less than 4.5 (soft drinks) and spoil the product • Microbiological results usually be obtained after 24 hours or at most two days 	204728547	600 g Alubag
TransFast-Gel 2 (TF-G 2 pH-Correction Fluid) Preparation: TF-G 1 + TF-G 2 Application: see above	204729782	9 x 250 ml
TransFast Gel Tubes (80 ml)	204730001	Bag (100 p.)

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NBB® media are containing heat-sensitive growth promoting substances for beer-spoiling bacteria. Therefore, a sparing, short as possible and not multiple heating resp. boiling should be applied.

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NBB® Culture Media for Beer-Spoiling Bacteria

of yeast sediments from bottles) are incubated in NBB-B (approx. 0.5 ml yeast suspension in approx. 10–20 ml NBB-B) and yeast-hazy beers (green beer, unfiltered beer, Bavarian or German wheat beer) are mixed with approx. 5 % NBB-C. For analysis of weak spots swabs are incubated in NBB-B-AM.

Due to the selectivity of beer all microorganisms growing clearly under correct conditions are relevant for biological production control. The harmless accompanying flora will be inhibited, especially hop-sensitive germs, e.g. lactic acid bacteria, which are not from the brewery. The selective parts of all NBB® media

incubation temperatures have to be between 25 and 28 °C and the final evaluation has to be made after 5 days for the NBB-A and the NBB-B samples, after 3 days for the NBB-B-AM samples and after 7–12 days for the NBB-C ones. Compared with other detection methods (several days) and normal shelf life samples (several weeks) you have a gain in time, especially with trace contaminations (1 germ in the sample) and difficult-to-detect beer-spoiling bacteria. You get an additional gain in time (1–2 days) by the improved composition of growth promoting substances and inhibitors and the process optimization in the new NBB® since 2001. As part of our continu-

The nutrient media NBB® has been specifically developed to detect beer-spoiling microorganisms. Depending on the types of samples collected in the brewery, a choice can be made among NBB-Agar (NBB-A), NBB-Broth (NBB-B), NBB-Broth-Enrichment-Medium (NBB-B-AM) and NBB-Concentrate (NBB-C). The different variations allow analysis of larger volumes and also of samples which are cloudy and/or contain large numbers of yeast cells: a detection of even minute traces of beer spoilers is thus possible. Rinse water and bright beer are for example incubated on NBB-A after membrane filtration.

Yeast samples (0.5–1.0 ml of pure culture, crop-pitching-tank-bottom yeast or

*are also suited to the sample type. The nutrients and the growth-promoting substances as well as the added stimulating compounds contained in NBB® are the features which give NBB® the capability for a quick and reliable detection of any beer spoiling bacteria using any procedure. As the quickness of the response depends on several factors (initial cell count, type of germs, degree of adaptation to beer, physiological condition and origin of the germs), an incubation as short as one day can be enough in case of higher contaminations, while mostly several days are needed in case of trace contaminations or of very slowly growing strains (e.g. *Lactobacillus lindneri*). Experience has shown that in order to attain a high degree of reliability the*

*ing efforts to refine and optimize our nutrient media, we have changed the inhibitory components of NBB-Broth and NBB-Concentrate. The change made to the liquid media was essential, because Actidione also slows the growth of *L. lindneri*, *Pediococcus damnosus* and other demanding sensitive beer-spoiling microorganisms. The pre-enrichment processes for rapid-identification methods by VIT and PCR also require a max. 2 days. With the earlier version, adequate pre-enrichment was not always possible or clear cut after an incubation period of up to 48 hours.*

*The major advantage of the new inhibitory components is that all beer-spoiling microorganisms (including *L. lindneri**

and *L. brevisimilis*) grow at equal rates (1–2 days time gain). In addition, sensitivity is enhanced, so that harmless accompanying GRAM-negative forms are much more effectively inhibited, making the findings more clear cut. Yeasts are likewise inhibited. The inhibitor dose was selected so that culture yeasts (as before) grow briefly and, as a result of CO₂-production, generate a typi-

cal anaerobic atmosphere in the enriched liquid. After a few hours they are then inhibited or killed. In the process, the yeast cells may take on atypical forms (involution forms). A disadvantage is the inhibition of GRAM-negative acetic acid bacteria and specific Enterobacteriaceae, since these play a very important role as indicator microorganisms in the analysis of weak spots, e.g. at the filling/capping machine or in the prefiltration area. For this reason, a suitable broth, called NBB-B-AM, was developed for swab samples. This provides optimum proliferation for all beer-spoiling microorganisms (including *Pectinatus* and *Megasphaera* spp.); moreover, it does not entirely inhibit the growth of acetic acid bacteria, other indicator microorganisms (e.g. *Enterobacter* spp.)

analysis of weak spots. With NBB-B-AM, the greater volume of CO₂ sometimes produced in the presence of yeasts may give rise to higher pressures in closed sample bottles. Care should therefore be taken to use burst-proof sample bottles (e.g. plastic sample bottles) in the analysis of weak spots. For the actual analysis of brewer's yeast, only normal NBB-B should be used.

media (NBB-B, NBB-B-AM) can be reduced. Usually 1–2 days are adequate, as the rapid-identification method works with a very low detection limit (100–1000 germs/ml). A reduction time for the pre-enrichment is also possible by VIT, based on the Fish-technology (fluorescence in situ hybridization). With specific marked gene probes and a fluorescence microscope for easy and reliable evaluation you receive


cal anaerobic atmosphere in the enriched liquid. After a few hours they are then inhibited or killed. In the process, the yeast cells may take on atypical forms (involution forms). A disadvantage is the inhibition of GRAM-negative acetic acid bacteria and specific Enterobacteriaceae, since these play a very important role as indicator microorganisms in the analysis of weak spots, e.g. at the filling/capping machine or in the prefiltration area. For this reason, a suitable broth, called NBB-B-AM, was developed for swab samples. This provides optimum proliferation for all beer-spoiling microorganisms (including *Pectinatus* and *Megasphaera* spp.); moreover, it does not entirely inhibit the growth of acetic acid bacteria, other indicator microorganisms (e.g. *Enterobacter* spp.)

NBB-B-AM is also suitable for analysis of weak spots for the beverage industry, as similar indicator microorganisms as in the brewery (acetic acid bacteria, *Klebsiella*, slime-shaping yeasts and moulds) are of particular relevance by producing biofilms. The new compositions of NBB-B and NBB-C do not alter the procedure or the predictive value of the tests in any way. Indeed, it has been found in practice that the selective detection of beer-spoiling microorganisms is even faster and yields more clear-cut results. Therefore, pre-enrichment as an application for rapid-identification methods by VIT (Vermicon Identification Technology) and PCR (Polymerase Chain Reaction) is confirmed.

Due to the following in-vitro specific amplification of the DNA by the PCR technique the pre-enrichment time in the selective

the results after max. 48 hours. Only the selected and living microorganisms will be specifically detected. These rapid-identification methods are based on molecular-biological fundamentals. A contamination can be explicitly detected and identified at the same time.

The easy-to-use bypass-membrane-system can be used for continuous sampling in the filtered area and for validation of filling machines for clear products. By obvious greater sampling volume the increase in detection security for trace contaminations is reached by using the common procedures for NBB-B, NBB-B-AM or NBB-C. For rapid-identification methods the bypass-membrane-system is the basis. Due to the described advantages a high security will be obtained.

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Combined Test Plan for a Secure Recording of Trace Contaminations

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Fig. 1: Diagramm of beer production and sampling

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Classification of Beer-Spoiling Bacteria

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Examination of Yeast Samples

*A very high biological purity is required for yeast samples. The goal is hence to detect not only absolutely harmful beer-spoiling organisms, but also all potentially harmful species, the indirectly beer-spoiling organisms and, in case of aerobic incubation, also the indicator germs. To the last mentioned belong in the first place the acetic acid bacteria (e. g. *Acetobacter pasteurianus* and *Gluconobacter frateurii*) as well as some GRAM-negative *Enterobacteriaceae* (e. g. *Enterobacter*-species, *Pantoea agglomerans*, *Obesumbacterium proteus*), Catalase-positive micrococci or lactococci. NBB-B is applied in such cases, as it's detection-sensitivity is the highest.*

If the contaminations are high enough and if the organisms are typical beer-spoilers, the colour of the indicator in the samples turns from red to yellow. In case of trace contaminations or of checks of older, autolyzing yeast cells from tank-bottoms, the change of the indicator colour is not evident enough, and the samples have to be examined


Abb. 4: Nachweis von bierschädlichen Bakterien in Hefeproben Fig. 4: Detection of beer-spoiling bacteria in samples of yeast

under the microscope before a final answer can be given. But as the growth of the beer spoiling bacteria is quite strong, the results are clear and easy to interpret, with several or even many suspect cells in each microscopic field. Scattered bacterial cells (e. g. 1-3 rods in 10 microscopic fields at a magnification of 600x) are however either dead or belong to harmless species, and were probably present in the (contaminated) yeast already before incubation. In doubtful cases it is advisable to transfer a sample, with a suitable inoculating loop, into new NBB-broth. If living beer-spoiling organisms are present, with this two-stage incubation they exhibit a strong growth and the colour of the indicator turns within 2 days.

The most suitable containers for these tests are test tubes with porous stoppers (e. g. plastic-foam plugs). In order not to stress unnecessarily the heat-sensitive culture medium it is recommended to use already sterilised containers, and to fill them, close to a flame, directly from the original bigger flasks. For larger samples of thick yeast slurry

(1-2 ml) also 50-ml-swing stopper sample bottles (50-ml-SSSBs) are quite suitable, into which 20-30 ml NBB-B are added (fig. 4). The examination of larger yeast samples (approx. 10 ml) could be done under the use of a shaker. A very practical way to check yeast samples are also sterile swabs (if possible with screw stopper). The sample of yeast or of tank bottoms is taken with the swab from various spots of the yeast container (total volume approx. 0.5 ml) and the tube is filled up to approx. 70 % with NBB-B. The screw-stopper should not be tight during the incubation, in order to allow the CO₂ to escape. If the samples are too large and/or if the yeast sediment is too dense, in case of trace contaminations there is the danger that some cells of the beer-spoiling organisms do not come in contact with the NBB-B. Under such circumstances it is advisable to re-suspend the sediment (e. g. by means of a slowly moving shaking table).

In checking yeast for purity anaerobic conditions are not absolutely necessary, as the yeast itself, at the start of the fermentation, provides for a favourable CO₂-environment in the sample. But after a short time the selective component present in the medium inhibits the culture yeast, so that the nutrients and the growth promoting substances are completely and without any restriction available for the beer spoilers, if they are present.

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Fig. 5: Selective detection of beer-spoiling bacteria corresponding to harmfulness. Add additional 5 % NBB-C of yeast containing samples. Filtered beer be alternatively admitted with approx. 10 % NBB-B.

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Examination of Samples of Unfiltered Beer

The purity of samples of green beer or of beer from the maturation stage or of wheat beer containing yeast can very easily be checked with NBB-concentrate (NBB-C). To the hazy samples, collected in 180-ml-swing stopper sample bottles (180-ml-SSSBs), 5 % (± 1 %) NBB-C is added and if necessary the bottles are filled up to the brim with sterile water. To avoid over-pressure the swing-stoppers should be closed tightly only later. These samples exhibit the highest beer-specific selectivity (comparable with shelf-life samples). But this selectivity can even be steered individually, by filling up the bottles more or less with sterile water (fig. 5). An anaerobic incubation is not necessary as the beer volume is high. Attention has, however, to be paid to use NBB-C always together with beer, as otherwise no beer-specific selectivity is given. It has often resulted advantageous to draw samples of approx. 120–150 ml green beer and to fill them up to the brim with sterile water after the addition of NBB-C. In this way, the higher concentrations of growth-inhibiting hop bitter substances are reduced. Furthermore, this procedure allows to obtain brim-full

bottles without any effort. A minimum of 7 and a maximum of 12 days at 27 °C are recommended for the incubation. While in these samples yeasts, and particularly culture yeasts, are inhibited (but NBB-C has to be added immediately to the sample containing yeast, or poured into the bottle before the sample is drawn), beer-spoiling bacteria grow very well in such a competition-free environment (fig. 6), and attain high germ counts, so that the microscopic examination becomes unproblematic (there is

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Fig. 6: Green beer sample
left: cloudiness by culture yeast
right: extra cloudiness by beer-spoiling microorganisms

no colour change in NBB-C). Also here, as we have already seen in discussing the evaluation of the NBB-B samples, only undoubted microscopic results should be considered: this requires that several suspect bacteria can be seen in each microscopic field. Scattered rods or cocci (e.g. 1 cell in 10 microscopic fields

at a magnification of 600x) are usually dead or harmless latent germs. In doubtful instances it is again advisable to carry out a second enrichment in NBB-B. Although in case of a contamination the hazy samples incubated with NBB-C exhibit usually an increase of the turbidity and/or of the sediments, their final microscopic evaluation is absolutely necessary. Some beer-spoiling organisms (e.g. *Pediococcus damnosus*) grow preferentially in the yeast sediment. In case of larger sample sizes (e.g. a “drop-flask” containing 1–2 liters of beer), it is advisable to proceed as follows: at first, the bright supernatant liquid is membrane-filtered, then the yeast sediment is suspended and poured into a SSSB into which the membrane filter and 5 % NBB-C are added before filling it up to the brim with sterile water.

Also the examination of unfiltered beer samples can be carried out simply and reliably with NBB-B. A large sample of unfiltered beer is drawn into a sterile container or into a suitable flask, which are kept until the yeast sediments to some extent (most of the beer-spoiling organisms which may be present sediment together with the yeast). The supernatant liquid is then poured away. With a Pasteur pipet or with a sterile swab approx. 0.5 ml yeast suspension are taken from various spots and the check is carried out as described under “Yeast Samples”. The selectivity is slightly reduced, but the time needed for a detection of the spoiling organisms is remarkably shorter.

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Examination of Filtered Beer Samples

Bright beer samples are usually membrane-filtered. The filter is incubated in Petri dishes on NBB-Agar at 27 °C, placed in anaerobic jars or in exsiccators blanketed with CO₂. For a better gas exchange the covers of the Petri dishes should be dented.

Beer-spoiling bacteria usually form plain colonies and turn the colour to yellow. Residues of oxygen in the membrane filters hinder however not seldom the growth of strictly anaerobic organisms like *Pectinatus* and *Megasphaera*, allowing at the same time a slight growth of the aerobic ones (mostly acetic acid bacteria).

Furthermore, air blisters in the membrane filters may impair the rapid diffusion of the nutrients and of the growth promoting substances. For these reasons the membrane filters are often incubated in NBB-B anaerobic jars. The culture flasks should be equipped with gas-permeable plugs (screw stoppers not to be tightened!).

Further ways to process the samples:

- Fill up to the standard filling level a sterile 50-ml-swing stopper sample bottle (50-ml-SSSB) with standard beer taken from bottles subjected to the biological shelf life test (without sterilizing it). Put the membrane filter into the sample bottle and fill up to the brim with NBB-B. This procedure allows an optimum detection of beer-spoiling organisms including harmful yeasts (haze, sediment, no colour

change). If the volume of the sample is large (drop flasks), pour approx. 40 ml into the sample bottle and membrane filter the rest, then put the membrane filter into the sample bottle and fill up to the brim with NBB-B.

- Take a beer bottle from the bottling line, fill it up to the brim adding 10–15 ml NBB-B and crown it again. Incubate it in the shelf-life cabinet for 10 days instead of 6 weeks. Also fermenting yeasts (both culture and foreign) can be detected. Haze, sediment, no colour change.
- Add to a sample of filtered beer 5 % NBB-C. Brim-full sample bottle (add sterile water) or anaerobic incubation. Detection of any beer-spoiling organisms, inhibition of yeasts.

Detection of Beer-Spoiling Bacteria in Samples of Wort and Water

50 ml of the sample (e. g. pitched wort or rinse water) are poured into a 180-ml-SSSB and 8 ml of NBB-C are added. In order to create a beer-specific, selective environment the sample bottle is then filled up to the brim with pasteurized beer and incubated for approx. 5 days at 27 °C.

Examination of Compressed Air and CO₂

At a sampling point the gas is allowed to slowly flow (1 min.) into NBB-B (e.g. 20 ml NBB-B in a SSSB). No preliminary flushing should be allowed, as otherwise also the beer-spoiling organisms which may be present are blown away.


Abb. 7: Wachstum von bierschädlichen Bakterien auf NBB-Agar c) *P. damnosus* d) *L. lindneri* (MF-Kultur Unterseite)
Fig. 7: Growth of beer-spoiling bacteria on NBB-Agar c) *P. damnosus* d) *L. lindneri* (MF-culture underside)

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The generation of weak spots in brewery and soft drinks industry is shown in fig. 11. Analyzing for indicators of beverage spoiling microorganisms will be prepared with NBB-swabs, analyzing for faecal indicators with Lactose-broth-swabs.

The routine analysis of weak spots is an essential part of an effective biological quality control program.

In the production area the weak spots of the beer circuit and contact areas (blind caps, pumps, CO₂-lines, by-passes and others) are checked with sterile swabs with approx. 10 ml NBB-B-AM (fig. 9) and aerobic incubation for a maximum of 3 days at 27 °C. Weak-spots are present if within the incubation time the indicator germs (mainly acetic acid bacteria) turn the colour to yellow. In an unfavourable case beer-spoiling organisms or faecal indicators may already be present in the mixture. Provided that detected contaminations are correctly eliminated, it is sufficient to check these areas every 3 or 6 months.

In the packaging area the weak spots are checked with sterile NBB-B-AM swabs in the direct and indirect contact area of the soaker discharge, of the bottle inspector and particularly of the filler and capper. Especially starwheels, conveyors, gears

and inner surfaces in the filler area are frequently niches of contamination and cause secondary infections (fig. 10). By experience it has been found that the hygienic conditions can be very well ascertained taking every week 20–30 swabs from each filling line (fig. 12).

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Special Samples (Swab)

1. Production/Beer Path

- Three-way cock
- CO₂-line
- Bunging apparatus
- CO₂-feeder
- Metering units
- Level indicator
- Blind caps
- Blind tubes, by-passes (especially on the separator, kieselgur filter and PVPP)
- Fittings
- Pumps
- Vessel sides (overflow channels)
- Vessel bottom
- Various installations and niches in the beer path
- Opened bottoms, sumps (hoses, tank outlets)
- Waste beer tank/old yeast tank (surroundings)
- Valve junction (underside)

Fig. 9: Special samples production area (swabs)

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Special Samples (Swabs)

2. Filling Area

Washer head (condensed water from different places)

Bottle Inspector

- Starwheels
- Bottle conveyor
- Other moist spots

- Starwheel, inside
- Baffles
- Panellings/girders
- Guideway of plastic rail

Filler

- Filling tube
- Control valves
- Bells
- Lifters
- Feed form
- Starwheel surface

Capper

- Guide rods
- Piston
- Crouner-inlet rail, plate
- Panelling
- Starwheels
- Bottle conveyor

Fig. 10: Special samples filling area (swabs)

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Fig. 12: Weak spot swabs on a filling line

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Fig. 13: Documentation of the weekly results in the bottling area

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If within a period of several weeks in average less than 30 % of the findings are positive with NBB-B-AM (yellow colouring after aerobic incubation for 3 days at 27 °C), the biological condition can be considered satisfactory. Higher numbers of positive results, or their trend to increase, point at stiff and persistent biofilms, in which sooner or later also true beer-harmful organisms or other problematic germs will start to adapt themselves and to grow (fig. 13).

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Analysis of Air-Borne Germs

The analysis of air-borne germs is a useful supplement of the swabs in the bottling area. For this type of checks the air-germ collector RCS of BIOTEST (63303-Dreieich/Germany) is very useful. Air samples are usually drawn in addition to the swab samples at the soaker-discharge, at the swab inspector and from the area around the filler/capper. In order to obtain clear results air suction times of approx. 2 minutes have to be used, in conjunction with NBB-air-germ-indicators (NBB-LKI). The aerobic incubation (3 days at 27 °C) detects mainly indicator germs, while typical beer spoilers are detected by the anaerobic one (5 days at 27 °C). Fields turned yellow are considered positive and expressed as a percentage of the total fields (34). After aerobic incubation the

yellow fields have in average to be less than 30 %. By experience it is then sure that a satisfactory level of protection against secondary contaminations has been attained. The results of the air-germ indicators are usually comparable with those of the swab samples. If they are conspicuously higher, this means that hidden weak spots are present.

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
Fig. 8: Air-germs indicators
 a) from washer, bottle discharge b) from filler c) from capper
 positive results of the fields 5–20 % (a,b,c: sanitary conditions o.k.)
 d) a lot of colonies with *L. brevis* (sanitary conditions critical)

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NBB[®] Media

Further development of the NBB[®] range: optimizing of the growth-promoting substances and inhibitors and of the detection of traces. Adaption on molecular-biological methods.

Product	Article-No.	Package
NBB[®]-Agar (NBB-A) Application: <ul style="list-style-type: none"> • Selective for beer-spoiling microorganisms. Inhibition of culture yeasts and harmless accompanying flora • Membrane filtration (beer and water samples) • Smears of beer-spoiling microorganisms • Yeast smears (yeast samples should better analyzed with NBB-B to detect beer-spoiling microorganisms) 	204709782	9 x 250 ml
NBB[®]-Broth (NBB-B) Application: <ul style="list-style-type: none"> • Selective for beer-spoiling microorganisms. Mostly inhibition of culture yeasts and harmless accompanying flora (incl. indicator germs) • Usual plant yeast samples, yeast sediments (0.5 ml sample for 10–20 ml NBB-B) • Sediment or separated product of unfiltered beer (e.g. green beer, wheat beer) • Incubation of membrane filters • For the quick pre-enrichment of traces of beer-spoiling bacteria for molecular-biological methods (VIT, PCR) • Add (5–20 %) to filtered beer. In completely filled bottles (e.g. swing stopper sample bottles) it is optimal for all beer-spoiling microorganisms incl. <i>Saccharomyces diastaticus</i>, <i>Pectinatus</i> and <i>Megasphaera</i> (culture yeast and indicator germs are suitable to a partly growth). Membrane filters can also be incubated in these bottles 	204710782	9 x 250 ml
NBB[®]-Broth-Enrichment-Media (NBB-B-AM) Application: <ul style="list-style-type: none"> • Optimal medium for all beer-spoiling microorganisms. Good growth of the indicator germs, especially acetic acid bacteria with aerobic incubation. Growth of yeasts (as well culture yeast, but slowly). Partly inhibition of harmless accompanying flora • Perfect for swabs to detect biofilms in the filling and production area (analysis of weak spots). Also for analysis of weak-spots in the whole beverage industry, because you can find similar indicator germs (especially acetic acid bacteria, <i>Klebsiella</i>, slime forming yeasts) growing by biofilm-formation • Incubation of membrane filters • These samples are very suitable for incubation of membrane filters (incl. bypass-membrane-system) and for the quick pre-enrichment traces of beer-spoiling bacteria for molecular-biological methods (VIT, PCR); also good detection of <i>Pectinatus</i> and <i>Megasphaera</i> • Add (5–20 %) to filtered beer (see NBB-B) 	204706782	9 x 250 ml
NBB[®]-Concentrate (NBB-C) Application: <ul style="list-style-type: none"> • Concentration for application 5 % (3–7 %) due to the whole volume of the sample. Only use diluted with beer samples or beer parts • Especially for detection of traces of beer-spoiling bacteria in unfiltered beer (yeast containing beer samples) (green beer, beer from the sampling cock, wheat beer). Inhibition of yeasts • Also as an additive for shelf-life samples and beer samples of the beer and filling path (in completely filled bottles and non diluted beer samples it is selective for obligatory beer-spoiling microorganisms). Inhibition of yeasts and accompanying flora • Rinse water and wort samples • Selectivity can adjust by diluting the beer sample with sterile water (see fig. 5) • Suitable for pre-enrichment for molecular-biological methods (VIT, PCR) 	204711782	9 x 250 ml
NBB[®]-B-Tubes Application: see NBB-B	204723646	box (20 p.)
NBB[®]-Air-Borne-Bacteria-Indicators (NBB-LKI) Application: <ul style="list-style-type: none"> • Using, additionally, the air-borne bacteria collector RGS from BIOTEST (63303 Dreieich/Germany), it is possible to detect beer-spoiling microorganisms and indicator germs in the most critical sectors of the bottling and kegging area, perhaps also in the production area 	204724244	box (50 p.)
NBB[®]-Powder (NBB-P) Application: see NBB-A and NBB-B	204716462	300 g Alubag
Sterile swabs	204725244	box (100 p.)

Focussed Product Portfolio

The Innovation Network.